Document No.: **470313372**  
 Revision: **A**  
 Change Order: **CO-035068**  
 Title: **PRIMARY LABEL INTERMEDIATE DISPENSER  
 CARTON TRANSFER-PACK CONTAINER WITH MALE  
 LUER ADAPTER, 1000 ML, K6R2247**

Page: 1 of 1  
 Last Modified: 16 Jul 2019 JML  
 Color: Black

**PROOFREADING INSPECTION / RELEASED ARTWORK**

Proofreading Approval:

Print Name \_\_\_\_\_ Signature \_\_\_\_\_ Date \_\_\_\_\_

**Adapter - 1000 mL**  
 (ko) **TRANSFER-PACK Container With Male Luer**

**Adapter - 1000 mL**  
 (en) **TRANSFER-PACK Container With Male Luer**

12x 포장단위: 12 unit / carton

(en) **TRANSFER-PACK Container With Male Luer Adapter - 1000 mL**

For the separation and storage of blood or blood components.

Store at Room Temperature. Protect from freezing. Avoid excessive heat.

Definition of “Room Temperature”:

“The temperature prevailing in a working area.”

Reference: United States Pharmacopeia, General Notices. United States Pharmacopeial Convention, Inc. 12601 Twinbrook Parkway, Rockville, MD.

Dispose of waste in appropriate biohazard container or according to local regulatory requirements.

(ko) **TRANSFER-PACK Container With Male Luer Adapter - 1000 mL**

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(en) **Instructions for fluid transfer: Use aseptic technique.**

1. Clamp a hemostat on the TRANSFER-PACK container tubing. If needed, clamp a second hemostat on tubing with female Luer connector.
2. Remove male Luer adapter cover from TRANSFER-PACK container and cover from appropriate female Luer connector.
3. Insert male Luer adapter into female Luer connector with a twisting motion until firmly seated.
4. Remove hemostat(s) to initiate fluid transfer.
5. When desired amount of fluid has been transferred, clamp hemostat(s) on TRANSFER-PACK container tubing and on tubing with female Luer connector.
6. Make three hermetic seals closely spaced, on TRANSFER-PACK container tubing as well as on the tubing with the female Luer connector. Carefully sever center seals using care to avoid fluid splatter. Remove hemostats.
7. Label container appropriately and store as required.

**Warning:**

**Container made of PVC with DEHP.** This container is of PVC material that used di-(2-ethylhexyl) phthalate (DEHP) as plasticizer, where DEHP is reported to affect the development of male sexual organs and sperm formation in young animal tests. Although infinitesimal amount of DEHP may be eluted in such PVC containers, there is very little or no risk from exposure to DEHP. Hence, as the risk that may appear from not using this medical device is much larger than the risk concerned due to DEHP, there is no need to avoid using it.

Phthalates are an integral part of the medical device formulation. Blood or blood components collected, prepared or transfused using product may contain low levels of DEHP and present certain risks. In general, the highest risk is associated with medical procedures involving infants, pregnant or nursing mothers, or requiring high volume to body weight blood infusion (such as trauma patients and neonate replacement transfusion). The use for infusion of blood or blood components with limited DEHP content must be weighed against the risks associated with not performing the medical procedure.

This product is designed for single use. Biological contamination of this product will result from normal use. Re-use of this product can expose patients or donors to this biological contamination. This exposure can result in adverse reactions, including severe illness and possibly death. This product is not designed for re-sterilization. Re-sterilization of this product may not completely eliminate biological contamination and may result in malfunction.

© 2019 Fresenius Kabi AG. All rights reserved.

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**Fenwal International, Inc.**