



EU Quality Management System Certificate (MDR)

Pursuant to Regulation (EU) 2017/745 on Medical Devices, Annex IX Chapters I and III
(Class I Devices in sterile condition, with measuring function or reusable surgical instruments)

No. G11 037258 0015 Rev. 00

Manufacturer: **Fresenius Kabi AG**
Else-Kröner-Str. 1
61352 Bad Homburg
GERMANY

SRN Manufacturer: DE-MF-000009273

The Certification Body of TÜV SÜD Product Service GmbH certifies that the manufacturer has established, documented and implemented a quality management system as described in Article 10 (9) of the Regulation (EU) 2017/745 on medical devices. Details on device categories covered by the quality management system are described on the following page(s).

The Report referenced below summarises the result of the assessment and includes reference to relevant CS, harmonized standards and test reports. The conformity assessment has been carried out according to Annex IX Chapter I and III of this regulation with a positive result. As applicable the involvement of the notified body is limited to the aspects relating to:

- establishing, securing and maintaining sterile conditions,
- conformity of the devices with the metrological requirements,
- reuse of the device, in particular cleaning, disinfection, sterilization, maintenance and functional testing and the related instructions for use.

The certified quality management system is subject to periodical surveillance by TÜV SÜD Product Service GmbH. All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with.

For details and certificate validity see: [www.tuvsud.com/ps-cert?q=cert:G11 037258 0015 Rev. 00](http://www.tuvsud.com/ps-cert?q=cert:G11_037258_0015_Rev_00)

Report No.: 713194775

Valid from: 2022-08-01

Valid until: 2027-07-31

Christoph Dicks
Head of Certification/Notified Body

Issue date: 2022-08-01



EU Quality Management System Certificate (MDR)

Pursuant to Regulation (EU) 2017/745 on Medical Devices, Annex IX Chapters I and III
 (Class I Devices in sterile condition, with measuring function or reusable surgical instruments)

No. G11 037258 0015 Rev. 00

Classification:	I
Device Group:	A020108 - ENTERAL FEEDING SYRINGES A0280 - SYRINGES - ACCESSORIES A030101 - INFUSION CONTROLLERS A030103 - ENTERAL FEEDING CONTROLLERS A040101 - ADMINISTRATION AND ASPIRATION FILTERS A060399 - FLUID COLLECTION BAGS AND SYSTEMS - OTHER A0706 - ANTISEPTIC CONNECTOR PROTECTORS A0799 - ADAPTERS, CONNECTORS, RAMPS, STOPCOCKS, CAPS - OTHER A080299 - INFUSION BAGS AND CONTAINERS, SINGLE-USE - OTHER G0280 - GASTROINTESTINAL TUBES - ACCESSORIES V020380 - NEWBORN NUTRITION DEVICES - ACCESSORIES V020399 - NEWBORN NUTRITION DEVICES - OTHER
Device Properties:	MDS 1005.1 - Ethylene Oxide sterilization

The validity of this certificate depends on conditions and/or is limited to the following: none